

SAM-TECH DIAGNOSTICS

LOC-200

Biochemistry analyzer

Overview

Fully automatic biochemistry analyzer using micro-fluidic reagent discs to deliver lab-quality results from just three drops of blood.

Product Description

LOC-200 Automatic Biochemistry Analyzer is designed to perform the blood biochemistry and electrolyte item analysis(35 parameters), delivering reliable and accurate results in 8-12 minutes by 3 drops of blood sample(90-120?L).

Key Features

- ? 8-12 min per test
- ? 100 µL sample size
- ? Micro-fluidic disc
- ? Intelligent IQC
- ? No centrifugation
- ? No calibration
- ? 7" touch screen
- ? Weighs 4.5 kg
- ? Multi-language GUI

- ? Built-in printer
- ? 35 parameters
- ? 8 panels including Liver (ALB, ALT, AST), Renal (CRE, UREA, UA), Electrolytes (K+, Na+, Cl-), Lipid (CHOL, TG, HDL-C), and Myocardial enzymes.

Product Information

Price: KES 0

For More Information

Contact our sales team for detailed specifications and pricing

Email: info@sam-techdiagnostics.com

Phone: +254 123 456 789

Website: www.sam-techdiagnostics.com

© 2026 SAM-Tech Diagnostics. All rights reserved.