

SAM-TECH DIAGNOSTICS

MQ 60 ProB

CLIA System Analyzers

Overview

High-performance automatic immunoassay system featuring an integrated centrifuge for rapid whole blood processing.

Product Description

The Hotgen MQ60 proB is a fully-automated bench-top chemiluminescence immunoassay (CLIA) system designed for high-volume clinical laboratories, offering a rapid High Throughput of 48 tests per hour, which is nearly double the capacity of the standard MQ60 smart model. This is achieved through its Dual Module Flexibility, featuring two independent detection modules with a total of 12 channels that support both efficient batch testing and urgent STAT modes. The analyzer is equipped with an Enhanced Storage capacity of over 10 million patient results for high-volume record keeping and features a large, upgraded 13.3-inch touchscreen display for clearer status monitoring. It also boasts a broad and comprehensive Highest test menu in the MQ60 series, supporting over 100 diagnostic parameters for a wide range of conditions, including inflammation, cardiac markers, and tumor markers

Key Features

- ? Built-in centrifuge
- ? Whole blood direct
- ? 15-inch touch screen
- ? 24 sample positions
- ? No hidden costs

- ? Result in 17 mins
- ? Single-test format
- ? Internal barcode
- ? Compact design
- ? Optional UPS
- ? Over 100 test menu

Product Information

Price: KES 0

For More Information

Contact our sales team for detailed specifications and pricing

Email: info@sam-techdiagnostics.com

Phone: +254 123 456 789

Website: www.sam-techdiagnostics.com

© 2026 SAM-Tech Diagnostics. All rights reserved.