

# SAM-TECH DIAGNOSTICS

## MQ 60 smart

CLIA System Analyzers

### Overview

---

Fully automated bench-top CLIA analyzer utilizing mono-test cartridges for flexible and precise immunoassay testing.

### Product Description

---

MQ60 series mono-test automatic chemiluminescence immunoassay analyzers adopt chemiluminescence enzyme immunoassay technology with magnetic particle technology, through detecting the analyte substrate luminous intensity, to calculate the analyte concentration.

### Key Features

---

- ? 15-inch touch screen
- ? 24 sample positions
- ? No liquid circuit
- ? Single cartridge test
- ? Result in 17 mins
- ? Whole blood testing
- ? Built-in calibration
- ? LIS/HIS connection
- ? Maintenance-free

? Integrated UPS

? Over 100 test menu

## Product Information

---

**Price:** KES 0

### For More Information

Contact our sales team for detailed specifications and pricing

**Email:** [info@sam-techdiagnostics.com](mailto:info@sam-techdiagnostics.com)

**Phone:** +254 123 456 789

**Website:** [www.sam-techdiagnostics.com](http://www.sam-techdiagnostics.com)

© 2026 SAM-Tech Diagnostics. All rights reserved.