

SAM-TECH DIAGNOSTICS

MQ60 Plus

CLIA System Analyzers

Overview

Next-generation small-molecule immunoassay analyzer using sandwich method testing for high-throughput laboratory efficiency.

Product Description

The Hotgen MQ60 Plus is a high-capacity, fully-automated bench-top chemiluminescence immunoassay (CLIA) system designed for central laboratories, featuring four independent detection units which enable the simultaneous processing of up to 24 tests at the same time. This flexibility allows for efficient batch testing with a high throughput. The analyzer uses single-dose, cartridge-based reagents to prevent contamination and eliminate the need for additional consumables like water or a drain connection. It supports a comprehensive menu of over 100 tests, including cardiac, inflammation, and tumor markers, and offers extensive data storage and LIS connectivity for seamless lab integration.

Key Features

- ? 150 tests per hour
- ? 30 sample positions
- ? 12 reagent spots
- ? First report: 17 min
- ? Clot detection
- ? 24h refrigeration

- ? 300 reaction cups
- ? Steel probe needles
- ? USB 3.0 interface
- ? LIS connectivity
- ? Over 100 test menu

Product Information

Price: KES 0

For More Information

Contact our sales team for detailed specifications and pricing

Email: info@sam-techdiagnostics.com

Phone: +254 123 456 789

Website: www.sam-techdiagnostics.com

© 2026 SAM-Tech Diagnostics. All rights reserved.