

SAM-TECH DIAGNOSTICS

Serbio Jr

Automated ESR analyzers

Overview

Automated benchtop ESR analyzer designed for direct determination of results from blood/EDTA samples.

Product Description

Gemini said SABIO JR Automated ESR Analyzer is a bench top system for small & midium size laboratories. It has complete integration in the hematology laboratories by the direct determination of ESR from blood/EDTA samples with its smart and innovative technology combining “ccd cam” and internationally accepted Westergren method. SABIO JR is completely cost effective, provides patient comfort and environment friendly since it is using same EDTA tubes coming from CBC tests. 16 sample loading positions, suitability to continuous loading for 19 minutes when you start running the instrument, being maintenance free, being environment-friendly, LCD touch screen and data transmission to LIS are just some of the advanced features of the analyzer which provides laboratory confort and productivity.

Key Features

- ? CCD camera technology
- ? No citrate tubes
- ? No reagents required
- ? 48 samples per hour
- ? First result: 20 min

- ? 16 loading positions
- ? 7" LCD touch screen
- ? Built-in printer
- ? LIS communication
- ? Maintenance-free

Product Information

SKU: RTK-BUNDLE

Price: KES 0

For More Information

Contact our sales team for detailed specifications and pricing

Email: info@sam-techdiagnostics.com

Phone: +254 123 456 789

Website: www.sam-techdiagnostics.com

© 2026 SAM-Tech Diagnostics. All rights reserved.